


VARILUX®
#1

PROGRESSIVE LENS BRAND
RECOMMENDED BY OPTOMETRISTS
AND OPTICIANS*

VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective MARCH 2022



* Survey conducted in 2018 by an independent third party. Sponsored by Essilor. Results were reported by independent opticians and optometrists.

VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective MARCH 2022

VARILUX® X DESIGN		VARILUX® X FIT											VARILUX® X 4D		
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	ADD POWER*	MAX CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00
THIN & LITE®	★ Clear	1.67											-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓	✓								-12.00 to +9.00	0.75 to 4.00D	-6.00
	★ Clear**	1.74											-14.00 to +13.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™**	1.74	✓	✓									-14.00 to +13.00	0.75 to 4.00D	-6.00

VARILUX® S DESIGN		VARILUX® S FIT											VARILUX® S 4D		
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	ADD POWER*	MAX CYLINDER
PLASTIC	Clear	1.50											-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	TREXA® (TRIVEX®)	Clear	1.53											-10.00 to +6.00	0.75 to 4.00D
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	0.75 to 4.00D	-6.00
	THIN & LITE®	Clear	1.67										-12.00 to +9.00	0.75 to 4.00D	-6.00
THIN & LITE®	Transitions® Signature® GEN 8™	1.67	✓	✓	✓								-12.00 to +9.00	0.75 to 4.00D	-6.00
	Clear**	1.74											-14.00 to +13.00	0.75 to 4.00D	-6.00
	Transitions® Signature® GEN 8™**	1.74	✓										-14.00 to +13.00	0.75 to 4.00D	-6.00

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

VARILUX® PHYSIO® W3+		TEMPORAL F	NASAL H	VARILUX® PHYSIO® W3+ FIT											TEMPORAL F	NASAL H	VARILUX® PHYSIO® W3+ EYECODE™				TEMPORAL E	NASAL H
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER						
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Transitions® XTRActive® new generation	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00						
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓	SM						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00						
	Transitions® XTRActive® new generation	1.59	✓	SM	✓	SM							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00						
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.60	✓	✓			✓						-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00						
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00						
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00						
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00						
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓	SM						-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00						
	Transitions® XTRActive® new generation	1.67	✓	SM	✓	SM	✓						-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00						
	Clear**	1.74											-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00						
Transitions® Signature® GEN 8™**	1.74	✓										-18.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00							

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

***THIN & LITE* 1.74 lenses come systematic with choice of Crizal® Preveniria®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

★ Available in **Essential Blue Series™** lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

VARILUX® COMFORT MAX		VARILUX® COMFORT MAX FIT															
		LENS ID												LENS ID			
		€												FE			
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER	
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.53	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00	
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00	
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-6.00	

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

VARILUX® LENSES AVAILABILITY for Eyecare Professionals

Effective MARCH 2022

VARILUX COMFORT® W2+		TEMPORAL F	NASAL H	VARILUX COMFORT® W2+ FIT										TEMPORAL F	NASAL H	VARILUX COMFORT® W2+ EYECODE™				TEMPORAL E	NASAL H
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER					
PLASTIC	★ Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
	Essilor Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
	Xperio UV™ Polarized	1.50	✓	✓			✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
	Transitions® XTRActive® new generation	1.50	✓										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00					
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00					
	Essilor Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00					
	Xperio UV™ Polarized	1.59	✓	✓			✓						-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00					
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00					
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00					
THIN & LITE®	★ Clear	1.60				✓							-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00					
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00					
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00					
	Essilor Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00					
	Xperio UV™ Polarized	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00					
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00					
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-6.00				

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

VARILUX PHYSIO DRx™

LENS ID



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓								-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® Polarized™	1.50	✓										-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75, 10.50	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® Polarized™	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50, 4.00D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® XTRActive® Polarized™	1.67	✓										-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
Transitions® Signature® GEN 8™**	1.74	✓										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00	

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

***For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.


★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)



VARILUX PHYSIO SHORT DRx™

LENS ID



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓								-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	0.75 to 3.00D	-6.00
	Transitions® XTRActive® Polarized™	1.50	✓										-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75, 10.50	0.75 to 3.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00
	Transitions® XTRActive® Polarized™	1.59	✓										-10.00 to +8.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00, 7.50, 8.50	0.75 to 3.00D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Transitions® XTRActive® Polarized™	1.67	✓										-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00
Transitions® Signature® GEN 8™ **	1.74	✓										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00	

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

***For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)

VARILUX COMFORT DRx™

LENS ID



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER	
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.50	✓	✓	✓								-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® Polarized™	1.50	✓										☼	-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75, 10.50	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.50D	-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® Polarized™	1.59	✓										☼	-10.00 to +8.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00, 7.50, 8.50	0.75 to 3.50D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.50D	-6.00	
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00	
	Transitions® XTRActive® Polarized™	1.67	✓										☼	-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	0.75 to 3.50D	-6.00
	★ Clear**	1.74												-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00
Transitions® Signature® GEN 8™**	1.74	✓											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50D	-6.00	

MFH=17mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

***THIN & LITE* 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

***For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.


★ Available in **Essential Blue Series™** lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ **Essilor Sun Colors (solid or gradient): gray mist, sienna brown, forest gray, grape, plum, emerald blue)**



VARILUX COMFORT SHORT DRx™



MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER	
PLASTIC	★ Clear	1.50											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00	
	Essilor Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00	
	Xperio UV™ Polarized	1.50	✓	✓			✓	✓	✓			☼	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	0.75 to 3.00D	-6.00	
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00	
	Transitions® XTRActive® new generation	1.50	✓	✓	✓								-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	0.75 to 3.00D	-6.00	
	Transitions® XTRActive® Polarized™	1.50	✓											-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75, 10.50	0.75 to 3.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53											-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00	
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00	
	Transitions® XTRActive® new generation	1.53	✓	✓	✓								-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	0.75 to 3.00D	-6.00	
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59											-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00	
	Essilor Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	0.75 to 3.00D	-6.00	
	Xperio UV™ Polarized	1.59	✓	✓			✓	✓	✓	✓	✓	☼	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00/7.50, 8.50	0.75 to 3.00D	-6.00	
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00	
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓								-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	0.75 to 3.00D	-6.00	
	Transitions® XTRActive® Polarized™	1.59	✓											-10.00 to +8.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00, 7.50, 8.50	0.75 to 3.00D	-6.00
THIN & LITE®	★ Clear	1.60											-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00	
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	0.75 to 3.00D	-6.00	
	★ Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00	
	Essilor Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00	
	Xperio UV™ Polarized	1.67	✓	✓					✓				-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	0.75 to 3.00D	-6.00	
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00	
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00	
	Transitions® XTRActive® Polarized™	1.67	✓											-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00	0.75 to 3.00D	-6.00
	★ Clear**	1.74											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00	
Transitions® Signature® GEN 8™ **	1.74	✓											-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.00D	-6.00	

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Preveincia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

***For the complete list of Transitions color availability please refer to the Transitions Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

★ Available in Essential Blue Series™ lenses.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

SM Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.

☼ Essilor Sun Colors (solid or gradient: gray mist, sienna brown, forest gray, grape, plum, emerald blue)

VARILUX PHYSIO ENHANCED™		VARILUX PHYSIO ENHANCED FIT™											VARILUX PHYSIO ENHANCED EYECODE™			
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
TREXA® (TRIVEX®)	Clear	1.53											-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.53	✓	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓	✓									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear	1.67											-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓	✓	✓							-12.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74											-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00
	Transitions® Signature® GEN 8™**	1.74	✓										-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Previncia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

VARILUX® PHYSIO®		LENS ID H → e											POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER	
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC						
PLASTIC	Clear	1.50												-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓										-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓										-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
TREXA® (TRIVEX®)	Clear	1.53												-10.00 to +6.00	78/83	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59												-10.00 to +6.00	80/85	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓										-10.00 to +6.00	75/80	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓										-10.00 to +6.00	75/80	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50D	-6.00
THIN & LITE®	Clear	1.60												-12.00 to +8.00	78/83 73/78	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50, 4.00D	-6.00
	Clear	1.67												-12.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.50, 4.00D	-6.00
	Clear**	1.74												-18.00 to +9.00	68 78/84 73/79	1.75, 2.75, 3.25, 5.25, 6.50 8.00	0.75 to 3.50, 4.00D	-6.00

MFH=17mm

VARILUX PHYSIO SHORT™		LENS ID → e S											POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER	
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC						
PLASTIC	Clear	1.50												-10.00 to +6.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
TREXA® (TRIVEX®)	Clear	1.53												-10.00 to +6.00	81/86	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59												-10.00 to +6.00	80/85	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.00D	-4.00
THIN & LITE®	Clear	1.60												-12.00 to +8.00		1.75, 3.00, 4.00, 5.50, 7.00 8.00	0.75 to 3.00D	-4.00
	Clear	1.67												-12.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50 8.00	0.75 to 3.00D	-4.00

MFH=14mm

NOTES: *ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ Currently Available

VARILUX COMFORT ENHANCED™		VARILUX COMFORT ENHANCED FIT™														
MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.60	✓										-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Clear	1.67											-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00
	Transitions® Signature® GEN 8™	1.67	✓	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00

MFH=14mm

NOTES: All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare lens treatment is applied (for example, Crizal® No-Glare lenses).

*ADD Power in steps of 0.25D.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

VARILUX COMFORT® 

MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +7.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +7.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-4.00
	Transitions® XTRActive® new generation	1.50	✓										-10.00 to +6.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-4.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Essilor Polarized	1.59	✓										-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Xperio UV™ Polarized	1.59	✓										-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
	Transitions® XTRActive® new generation	1.59	✓										-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-4.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.60	✓										-12.00 to +8.00	73/78	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00
	Clear	1.67											-14.00 to +9.00	78/83	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-4.00
	Transitions® Signature® GEN 8™	1.67	✓	✓									-14.00 to +9.00	73/78	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00

MFH=17mm

VARILUX COMFORT SHORT™ 

MATERIAL	PRODUCT	INDEX	G	B	GG	COLORS	G1	B1	G-G	C	R	ESC	POWER RANGE	DIAMETER	BASE CURVES	ADD POWER*	CYLINDER
PLASTIC	Clear	1.50											-10.00 to +7.00	80/85	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.50	✓	✓									-10.00 to +7.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
	Transitions® XTRActive® new generation	1.50	✓										-10.00 to +6.00	77/82	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.00D	-4.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59											-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Essilor Polarized	1.59	✓	✓									-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Xperio UV™ Polarized	1.59	✓	✓									-10.00 to +6.00	81/86	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.59	✓	✓	✓								-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
	Transitions® XTRActive® new generation	1.59	✓										-10.00 to +6.00	78/83	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.00D	-4.00
THIN & LITE®	Clear	1.60											-12.00 to +8.00	78/83 73/78	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.60	✓										-12.00 to +8.00	78/83 73/78	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.00D	-4.00
	Clear	1.67											-14.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-4.00
	Transitions® Signature® GEN 8™	1.67	✓	✓									-14.00 to +9.00	78/83 73/78	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.00D	-4.00


MFH=14mm

NOTES: *ADD Power in steps of 0.25D.

G: Gray, B: Brown, GG: Graphite Green, Colors: Transitions® Signature® style colors in Sapphire, Amethyst, Amber, or Emerald, G1: Gray 1, B1: Brown 1, G-G: Gray-Green, C: Copper, R: Ruby, ESC: Essilor Sun Colors

✓ **Currently Available**

VARILUX® LENS PORTFOLIO



VARILUX® LENS TECHNOLOGIES

PATIENT BENEFITS



©2022 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries. *Transitions, Transitions Signature,* and *XTRActive* are registered trademarks and the *Transitions* logo is a trademark of Transitions Optical, Inc., used under license by Transitions Optical Ltd. GEN 8 is a trademark of Transitions Optical Limited. **Photochromic performance is influenced by temperature, UV exposure, and lens material.** These products may be protected by one or more patents listed at www.essilorusa.com/patents. 209954_PRO_VAR_LL/ECST 03/22